



Green button
75 mg




Grey button
150 mg

Label
Colour coded according to strength and includes use by date

Injection button
(green or grey)
Audible clicks at beginning and end of injection²

Window
To visualise drug delivery²

Yellow safety cover
Needle not visible to patient
Latex free²



Blue cap

All patients agreed that learning to use the PRALUENT pre-filled pen is easy because all the parts are clearly labelled and colour-coded**⁵

- ✓ Self-injected or caregiver administered¹
- ✓ Same device available in two dosing strengths (75 mg and 150 mg)¹